



Recombinant Probable glutathione S-transferase 5 (gst-5)

Product Code	CSB-BP600253CXY
Abbreviation	gst-5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09596
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVSYKLTYFN GRGAGEVSRQ IFAYAGQQYE DNRVTQE QWP ALKETCAAPF GQLPFLEVDG KKLAQSHAIA RFLAREFKLN GKTAWEEAQV NSLADQYKDY SSEARPYFYA VMGFPGDVE TLKKDIFLPA FEKFGFLVN FLKASGSGFL VGDSLWIDL AIAQHSADLI AKGGDFSKFP ELKAHAEKIQ AIPQIKKWIE TRPVTPF
Source	Baculovirus
Target Names	gst-5
Protein Names	Recommended name: Probable glutathione S-transferase 5 EC= 2.5.1.18 Alternative name(s): GST class-sigma
Expression Region	1-207
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.