



Recombinant *Oryza sativa* subsp. *indica* Abscisic acid 8'-hydroxylase 2 (CYP707A6), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP600313OFF |
| Abbreviation | CYP707A6 |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q09J78 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. <i>indica</i> (Rice) |
| Purity | >85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | CYP707A6 |
| Protein Names | Recommended name: Abscisic acid 8'-hydroxylase 2 Short name= ABA 8'-hydroxylase 2 EC= 1.14.13.93 Alternative name(s): Cytochrome P450 707A6 OsABA8ox2 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |