



Recombinant *Saccharomyces cerevisiae* Delta (3,5)-Delta (2,4)-dienoyl-CoA isomerase (DCI1)

Product Code	CSB-YP600503SVG
Abbreviation	DCI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08558
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSSRVCYHIN GPFFIIKLID PKHLNSLTFE DFVYIALLLH KANDIDSVLF TVLQSSGKYF SSGGKFSAVN KLNDGDVTSE VEKVSCLVSA ISSPNIFVAN AFAIHKKVLV CCLNGPAIGL SASLVALCDI VYSQNDSVFL LFPFSNLGFV AEVGTSVTLT QKLGINSANE HMIFSTPVLV KELIGTIITK NYQLTNTETF NEKVLQDIKQ NLEGLYPKSV LGMKELLHSE MKQKLIKAQA METNGTLPFW ASGEPFKRFK QLQEGNRRHK L
Source	Yeast
Target Names	DCI1
Protein Names	Recommended name: Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase EC= 5.3.3.-
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.