



Recombinant *Enterobacter aerogenes* HTH-type transcriptional regulator cbl (cbl)

Product Code	CSB-EP600505EKN
Abbreviation	cbl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08598
Product Type	Recombinant Protein
Immunogen Species	<i>Klebsiella aerogenes</i> (<i>Enterobacter aerogenes</i>)
Purity	>85% (SDS-PAGE)
Sequence	MNFQQLKIIR EAARQDYNLT EVANMLYTSQ SGVSRHIREL EEELGIEIFI RRGKRLLGMT EPGKALLSIA ERILNEASNV RRLADLFTND ASGVLTIATT HTQARYSLPP VIKAFRELFS DVRVELVQGT PQEIEALLHN GGADIGIASE RLSNDPTLAA FPWFRWHHSL LVPKDHPLTQ VSPLTLEAIA RWPLITYRQG ITGRSRIDEA FNRKGLMPDI VLSAQSDSVI KTYVELGLGV GLVAEQSGDA READTFTRLD TRHLFDANTV WLGLKRGQLQ RNYVWRFIEL CNAGLSLDEI KRQAMEPEEV AIDYQI
Source	<i>E.coli</i>
Target Names	cbl
Protein Names	Recommended name: HTH-type transcriptional regulator cbl
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.