



# Recombinant Guinea pig C-C motif chemokine 2 (CCL2)

<b>Product Code</b>	CSB-MP600526GU
<b>Abbreviation</b>	CCL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q08782
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Cavia porcellus (Guinea pig)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QPDGVNT PTCCYTFNKQ IPLKRVKGYE RITSSRCPQE AVIFRTLKNK EVCADPTQKW VQDYIAKLDQ RTQQKQNSTA PQTSKPLNIR FTTQDPKNRS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CCL2
<b>Protein Names</b>	Recommended name: C-C motif chemokine 2 Alternative name(s): Monocyte chemoattractant protein 1 Monocyte chemotactic protein 1 Short name= MCP-1 Small-inducible cytokine A2
<b>Expression Region</b>	24-120
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene is one of several cytokine genes clustered on the q-arm of chromosome 17. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. This protein is structurally related to the CXC subfamily of cytokines. Members of this subfamily are characterized by two cysteines separated by a single amino acid. This cytokine displays chemotactic activity for monocytes and basophils but not for neutrophils or eosinophils. It has been implicated in the pathogenesis of diseases characterized by monocytic infiltrates, like psoriasis, rheumatoid arthritis and atherosclerosis. It binds to chemokine receptors CCR2 and CCR4.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.