



Recombinant Danio rerio Regulator of G-protein signaling 7-binding protein A (rgs7bpa)

Product Code	CSB-MP600565DIL
Abbreviation	rgs7bpa
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08BU8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSSAPNGRKN RPRTAGTIFQ IGGKAPSRES ERRESSEALL AVSDCRMVVQ EFNTLVVALYR ELVISIGEIS ADCPSLRAEM HKTRTKGCEM ARTAHQNLSA ISGPEDGEIH PEICRLFQIL QCCELYLMLTE MLKSVCLLGS LQLHRKGKHY LGATKLDNRK EDSSDIPILE DTSSTPPDCP HTYFLVATDI ENIERDMTEM KNLLSKLRET MPLPLKNQDD SLLNLAPYP MVRQRKRRFF GLCCLVSS
Source	Mammalian cell
Target Names	rgs7bpa
Protein Names	Recommended name: Regulator of G-protein signaling 7-binding protein A Alternative name(s): R7 family-binding protein A
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.