



# Recombinant Bovine Endophilin-B2 (SH3GLB2)

<b>Product Code</b>	CSB-EP600575BO-B
<b>Abbreviation</b>	SH3GLB2
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q08DK5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MDFNMKKLAS DAGIFFTRAV QFTEEKFGQA EKTELDAHFE SLLARADSTK NWTEKILRQT EVLLQPNPSA RVEEFLYEKL DRKVPSRVTN GELLAQYMAE AASELGPTTP YGKTLIKVAE AEKHLGAAER DFIHTASINF LTPLRNFLEG DWKTISKERR LLQNRRLDLD ASKARLKKAK AAEAKATTVP DFQETRPRNY ILSASASALW NDEVDKAEQE LRVAQTEFDR QAEVTRLLLE GISSTHVNHL RCLHEFIESQ TTYAQCYRH MLDLQKQLGR FPGTFVGTAE PASPPLSSTS PTTTAATMPM GPSVADLAPP GEAALRLEEV APPASGTRKA RVLYDYEAAD SSELALLADE LITVYSLPGM DPDWLIGERG NKKGKVPVTV LELLS</p>
<b>Source</b>	E.coli
<b>Target Names</b>	SH3GLB2
<b>Protein Names</b>	Recommended name: Endophilin-B2 Alternative name(s): SH3 domain-containing GRB2-like protein B2
<b>Expression Region</b>	1-395
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>