



# Recombinant Pig 5'-AMP-activated protein kinase subunit gamma-1 (PRKAG1)

<b>Product Code</b>	CSB-BP600597PI
<b>Abbreviation</b>	PRKAG1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q09138
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METVTSSDSS SAVENEHPQD TPESNNSVYT SFMKSHRCYD LIPTSSKLVV FDTSLQVKKA FFALVTNGVR AAPLWDSKKQ SFVGMLTITD FINILHRYK SALVQIYELE EHKIETWREV YLQDSFKPLV CISPNASLFD AVSSLIRNKI HRLPVIDPES GNTLYILTHK RILKFLKFI TEFKPEFMS KSLEELQIGT YANIAMVRTT TPVYVALGIF VQHRVSALPV VDEKGRVVDI YSKFDVINLA AEKTYNNLDV SVTKALQHRS HYFEGVLKCY LHETLETIIN RLVEAEVHRL VVVDENDVVK GIVSLSDILQ ALVLTGGEKP
<b>Source</b>	Baculovirus
<b>Target Names</b>	PRKAG1
<b>Protein Names</b>	Recommended name: 5'-AMP-activated protein kinase subunit gamma-1 Short name= AMPK gamma1 Short name= AMPK subunit gamma-1 Short name= AMPKg Alternative name(s): 38 kDa subunit
<b>Expression Region</b>	1-330
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.