



Recombinant Serpentine receptor class alpha-3 (sra-3)

Product Code	CSB-BP600604CXY
Abbreviation	sra-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09205
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSSQKCASHLEIARLES LNFKISQLIYFVLIITLFFTTFFALKVVIQKKCIFQLSTKILL Y QNLFSANIHQIFLGITIVERLNI AFFKLNDKCIPLRPETNCQLYLEMFIAGLSGMV YGQT GLLFERACATFIKKYEKRKSLATGIITSIIVLLLSGSTARIIIWDDPLSEYQLACTS WPS KSRDRSTMFFSICTFISLNLVISLLIRRHNEKLEYSTPFVVGPRFRKREVIDSTS TICF LFFFQFIFMFIYSFGVLLGTIREIIGYEQYYFWVWVYTIPFIAASFPILLIYRIRY SN TNRIMIIKQFTSTKQTQEDHIKQMKNVWN
Source	Baculovirus
Target Names	sra-3
Protein Names	Recommended name: Serpentine receptor class alpha-3 Short name= Protein sra-3
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.