



Recombinant Probable histone deacetylase complex subunit SAP18 (C16C10.4)

Product Code	CSB-EP600608CXY
Abbreviation	C16C10.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09250
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSNGQLVSQV TVGHDKPLDR EKVCPMLLRV FVANNRHNPM SEYNSRNGGS VPPSELQMHT WMDCSLRELT SLIKEVNPDA RRKGTTFDFA IVQADRGSPR YILRDVGNTT NGERGIDDNK TLQQCKFEVG DFVDVAISLP SQGRRFNNRE QGDRFDHRQR QRSPIR
Source	E.coli
Target Names	C16C10.4
Protein Names	Recommended name: Probable histone deacetylase complex subunit SAP18 Alternative name(s): 18 kDa Sin3-associated polypeptide
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.