



Recombinant Uncharacterized protein F10B5.2 (F10B5.2)

Product Code	CSB-EP600613CXY-B
Abbreviation	F10B5.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09305
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGGALPPEIV DMYRNGAFL LFLGFPQASE FGIDYKSWKT GEKFMGLKMI PPGVHFVYCS IKSAPRIGFF HNFKAGEILV KKWNTSESETF EDEEVPTDQI SEKKRQLKNM DSSLAPYPYE NYRSWYGLTD FITADTVERI HPILGRITSQ AELVSLETEF MENAEKEHKD SHFRNRVDRE NPVRTRFTDQ HGLPIMKIRE GYEIRFQDIP PLTVSQNRVG IEYSDRLYRL LRALGGDWKQ LLAEMQIAFV CFLQGQVFEG FEQWKRIIHL MSCCPNSLGS EKELFMSFIR VLFFQLKECP TDFFDIVSR DNFLTTLTSM LFANVRDSAHA AGDDLKKTAK QFKQYLTNQF KWDFNCD
Source	E.coli
Target Names	F10B5.2
Protein Names	Recommended name: Uncharacterized protein F10B5.2
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.