



Recombinant Uncharacterized protein F58F12.3 (F58F12.3)

Product Code	CSB-YP600623CXY
Abbreviation	F58F12.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09400
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MHSQFGKHKSDISLLTPAPPSPKAFGLLYAAHWLGSTDPEIWERLETLANENQ CLIDLRD EKVPGTYETQLFVIGYAGNIVTVLKQLQDYIDEPNNFDRFDIDRRPRKETLKMC VPTDMF LSKFGVGLVAEGRITHYEQNRGVTIKTYAISQSQHVISVYSNENNHQLPLHDLV EYFAYP SQCYEKSFAQKIPYLVDNKSGLLNRFKDISKLDVIKTPRDQLWLMALDGMKTV DARESIL QALRENPDSASLYVRACKLEHLPSCKMEWLLHALQRDTLDDKMKIIKHAETYF RGPAILQ ILESAVLSLSNDIHIVRELAIAHEKSGNFEKSMQFLKDCIDNNSWNSTELTASMI HCEHA GANRTSKVYLEAMYSIFNFTTSEEAAALIEIVKVLVETDSPLSAKSLSLELMTRYA GTPTF EKTRQLFHSLFYTTCSKSDSERYDDTFIKHMAKWNKKQQEQK
Source	Yeast
Target Names	F58F12.3
Protein Names	Recommended name: Uncharacterized protein F58F12.3
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.