



Recombinant Uncharacterized zinc finger protein C28H8.9 (C28H8.9/C28H8.10)

Product Code	CSB-EP600631CXY-B
Abbreviation	C28H8.9/C28H8.10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09477
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MISENGYMDL MRNCIKWNSR QMEDRKKRLR FPYYEHQTAT AQRESRFTTR NVDHMYPNSD PNTVVQFATE RWKKSQSQPP SDAVEMSMFL RENPSIQVAI DHLTPQVVGP TTESVSDSSN DSTTIRPSRQ TQIKEEYRDD YVLDDLESPD EFGSDEDDWS SRKRRKGNLG PVQKATSSRK KVPTTRSSVS RLTPSRISIVK ETKYEEPEEK TYPCKDKCSAK YKSLAGLSYH QSYLHDQKSS QPLVKLLSPS IEISTSCDFC SGTAFMNKNT KLPEDLVSCH DCGRSGHPSC LNFNQNVTKI IKRSGWQCLE CKSCTICGTS ENDDKLLFCD DCDRGYHLYC LTPALEKAPD DEYSCRLCQV EFGDKASAPA KK
Source	E.coli
Target Names	dpff-1
Protein Names	Recommended name: Uncharacterized zinc finger protein C28H8.9
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.