



Recombinant Uncharacterized GTP-binding protein E02H1.2 (E02H1.2)

Product Code	CSB-MP600634CXY
Abbreviation	E02H1.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09523
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MIKRLNFSC RLFSTASTST SYEKTPIQPA TKCLQLAVIG APNVGKSLLT NSLIRCPLSA VSSKMDTTTR NISASICSDS TQLVFVDSPG AVSTSHVRQT MKKTSATSGD RVLQDPERAL QRAQHVLVVQ DSTAPGAYIH HRVLHMLHRY SHVPSILVMN KIDLVMRRSD LLPLVEILTN GQLSDNQQIS TKPAQIGRLG KSLSTNIQSS SSKFPSDEKW QSQFRELIQK PTWKCSYSET RSLFRTICGW SGFERVFFVS SLNGEGIDEL RDHLMSISPQ GEWKMQDGMP TGESAQQLCI DSIRAAVLDT TPSDVAYTVQ IRISEWEEQG EVLQIVGEIR CQKPRDGSLI IGKGGKRISE IGRRVNEHLH SLFQRQLYAR LIVTHNGKLI TQSK
Source	Mammalian cell
Target Names	E02H1.2
Protein Names	Recommended name: Uncharacterized GTP-binding protein E02H1.2
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.