



Recombinant *Schizosaccharomyces pombe* Probable ubiquitin carboxyl-terminal hydrolase 8 (ubp8)

Product Code	CSB-EP600654SXV
Abbreviation	ubp8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09738
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPGDVEGCQHLKLPADVENYQKICTQIFSCHFVPRRCSTCKRINKRSIRCLS CHSVGCL WGHGEEHAMEHTHMIGVDVKNGHTYCFGCQDYVYQTELETLRFKIKNIKAW QSDHKRLP EKYNQMVCLEAYRKYPPVCATAGLRGIQNLGATCFMSVILQSILHNPLVRNLFF SGFHTS TDCKRPTCMTCAIDDMFSSIYNSKNKSTFYGPTAVLNLMWKLSKSLCGYSQQ DGHEFFVY LLDQMHTESGGGTSMPTCPIHRIFSGSLKNVVTCLDCKKERVAVDPLMDISL DINEPTL QGCLERFVSKEKVQYSCHSCGSKNAIKQLVFDKLPPTICMQLKRFEQNNFAMS TKIDKQV SYPAFLRMRYNFNQDDVDYQLYSVVCHKGTLDTGHYIAYTYQYQNQWFLDDT TIVEVKES EVLNSQAYLLFYHERQILYSDEMTVKTEN
Source	<i>E.coli</i>
Target Names	ubp8
Protein Names	Recommended name: Probable ubiquitin carboxyl-terminal hydrolase 8 EC=3.4.19.12 Alternative name(s): Deubiquitinating enzyme 8 Ubiquitin thioesterase 8 Ubiquitin-specific-processing protease 8
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.