



Recombinant Schizosaccharomyces pombe tRNA wybutosine-synthesizing protein 2 (trm12)

Product Code	CSB-BP600666SXV
Abbreviation	trm12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09832
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MISFLVQKSE AKYWKDVVQN KNKFLKNIGI ISGPKLVDTY SNQFEADFLK KCVLIPTSYN KEDLPKDILD CPRIYLIEEE VKSVQTDEKL QIKVLRYSQ ILHTTVSEEQ YPVPSRFVIY PPLALLSLNS FGTKEWHLFF QSYANVDVKA GLFDIFAREW NVSHIAINEP IPVGDVMRRP LSLKPLYGDF GVLIDGNPSE QDFQKAFWVS CRQNGIYQTW APLYTMFSRG NSIEKARVLN FSYIKDEIIA DLYAGIGYFT FSYVKAGAST VFCWEINPWS VEALRRAALK NGWSIYVVRN GENYEFKVG T HRIVVFLESN VYAAERFSKM NISARHINLG MLP SSEK SWS TATSILKRES HSFIVHENV KDEDIETYSS EVNSCF SKML AKNTVCVTNC VKSFSPRVSH IVYDIETV
Source	Baculovirus
Target Names	trm12
Protein Names	Recommended name: tRNA wybutosine-synthesizing protein 2 Short name= tRNA-yW-synthesizing protein 2 Alternative name(s): Alpha-amino-alpha-carboxypropyl transferase TYW2 EC= 2.5.1.-
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.