



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C16C9.01c (SPAC16C9.01c, SPAC4G8.14c)

Product Code	CSB-YP600667SXV
Abbreviation	SPAC16C9.01c, SPAC4G8.14c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09839
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDYPSKDVEF VTTGMFIIDE IEFDGGKKIK DIIGGAGSFA LIGARLWTSE PESQQLGMIV DKGSDFPQVI LDELNDLHTG IVYRETPWRK TTYGYNRLRD NGFRLFYILT PKKPIRAPDL LETGLIYAKT IHIITNPKTF VSVCEELAKY GNLKNRPKIV WEPTPESCLP ENWHILQKAL KYCDIFSPNE VDSANLLGID ISEKSRAKE FVSAFQEFQI GHDSQGWWVL RNGSDGCLIG ACTSSTNATS DAVVPVREIL HLPVQMKNG SVLDTTGAGN AFCGAAILEY YRTGNIIIEAC AKATVAASFV IQQHGLPAHT KNSEGIDLWN GESVFQRLKE YYEFLASHEA SAIKELTKRG V
Source	Yeast
Target Names	SPAC16C9.01c
Protein Names	Recommended name: Uncharacterized protein C16C9.01c
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.