



Recombinant Schizosaccharomyces pombe

Replication factor C subunit 2 (rfc2)

Product Code	CSB-EP600668SXV
Abbreviation	rfc2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09843
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSFFAPRNKK TEQEAKKSIP WVELYRPKTL DQVSSQESTV QVLKKTLLSN NLPHMLFYGS PGTGKTSTIL ALSRELFGPQ LMKSRVLELN ASDERGISII REKVKSFAKT TVTNKVDGYP CPPFKIIILD EADSMTQDAQ AALRRTMESY ARITRFCLIC NYMTRIIDPL SSRCSKYRFK PLDNENMVKR LEFIAADQAV SMEPGVVNAL VECSSGDMRK AITFLQSAAN LHQGTPITIS SVEELAGAVP YNIIRSLDT AYTKNVSNIIE TLSRDVAAEG YSTGIILSQL HDVLLKEETL SSPVKYKIFM KLSEVDKRLN DGADETLQLL DLLSSISVVC
Source	E.coli
Target Names	rfc2
Protein Names	Recommended name: Replication factor C subunit 2 Short name= Replication factor C2 Alternative name(s): Activator 1 41 kDa subunit
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.