



# Recombinant Human Serine/threonine-protein phosphatase 2B catalytic subunit alpha isoform (PPP3CA)

<b>Product Code</b>	CSB-EP600840HU-B
<b>Abbreviation</b>	PPP3CA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q08209
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SEPKAIDPK LSTTDRVVKA VPFPPSHRLT AKEVFDNDGK PRVDILKAHL MKEGRLEESV ALRIITEGAS ILRQEKNLLD IDAPVTVCGD IHGQFFDLMK LFEVGGSPAN TRYLFLGDYV DRGYFSIECV LYLWALKILY PKTLFLLRGN HECRHLTEYF TFKQECKIKY SERVYDACMD AFDCLPLAAL MNQQFLCVHG GLSPEINTLD DIRKLDRFKE PPAYGPMCDI LWSDPLEDFG NEKTQEHFTH NTVRGCSYFY SYPAVCEFLQ HNNLLSILRA HEAQDAGYRM YRKSQTTGFP SLITIFSAPN YLDVYNNKAA VLKYENNVMM IRQFNCSPHP YWLPNFMDFV TWSLPFVGEK VTEMLVNVLN ICSDDELGSE EDGFDGATAA ARKEVIRNKI RAIGKMARVF SVLREESESV LTLKGLTPTG MLPSGVLSSG KQTLQSATVE AIEADEAIKG FSPQHKITSF EEAKGLDRIN ERMPPRRDAM PSDANLNSIN KALTSETNGT DSNGSNSSNI Q
<b>Source</b>	E.coli
<b>Target Names</b>	PPP3CA
<b>Protein Names</b>	Recommended name: Serine/threonine-protein phosphatase 2B catalytic subunit alpha isoform EC= 3.1.3.16 Alternative name(s): CAM-PRP catalytic subunit Calmodulin-dependent calcineurin A subunit alpha isoform
<b>Expression Region</b>	2-521
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.