



Recombinant Schizosaccharomyces pombe 60S ribosomal protein L10-A (rpl10a)

Product Code	CSB-BP600969SXV
Abbreviation	rpl10a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09127
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MARRPARCYR YCKNKYPYKPKS RYNRAVPDSK IRIFDLGRKR AGVDEFPLCI HLVSNEYEQI TSEALEAARI CANKYLVKIG GKDSFHLRVR AHPFHVVRIN KMLSCAGADR LQTGMRHAFG KPNGLVARVN IGQILMSVRT KDSSRATAIE ALRRCQYKFP GQQRIIVSKK WGFSQYARDE YIEKRSRGEI IPDGCYAKFL NKRGSLSQEKL DLFPEASFNL A
Source	Baculovirus
Target Names	rpl1001
Protein Names	Recommended name: 60S ribosomal protein L10-A Alternative name(s): QM protein homolog SpQM
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.