



Recombinant Schizosaccharomyces pombe Protein phosphatase 2C homolog 2 (ptc2)

Product Code	CSB-EP600973SXV
Abbreviation	ptc2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09172
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGQTLSEPVL DKHSSSGGDR WLHFGVSHMQ GWRISMEDAH CALLNFTDSN SSNPPTSFFG VFDGHGGDRV AKYCRQHLPD IIKSQPSFWK GNYDEALKSG FLAADNALMQ DRDMQEDPSG CTATTALIVD HQVIYCANAG DSRTLVLGRKG TAEPLSFDHK PNNDVEKARI TAAGGFIDFG RVNGSLALSR AIGDFEYKKD SSLPPEKQIV TAFPDVVIHN IDPDDEFILIL ACDGIWDCKS SQQVVEFVRR GIVARQSLEV ICENLMDRCI ASNSESCGIG CDNMTICIVA FLHGRGLEDW YNWITQRVNS GEGPCAPPSY AELRGPNTIA DARNLQLEYD HIASHEYGSG DTYDSDSDDE TIAYDRYYLH
Source	E.coli
Target Names	ptc2
Protein Names	Recommended name: Protein phosphatase 2C homolog 2 Short name= PP2C-2 EC= 3.1.3.16
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.