



Recombinant Uncharacterized WD repeat-containing protein F21H12.1 (F21H12.1)

Product Code	CSB-EP600989CXY
Abbreviation	F21H12.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09309
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MVEILDRTYGAQFPPEELECHLDLQNASANCKFNRWGSIVAVGCTDGRVLIYD FMTRNIA RTFSAHCLPVSCLSWSRDGRKLLTSSADNSIAMFDVLAGTLLHRIRFNSMVTF AMFHPRN DNKAIVLQVKNQPTVEQFSPRIQTVLANDTPGSSDESASCVSYDRKGGKYIIAGT GKGKLL IYNAETLKCVAWCKQNTVQQIRQIIVPMKSRFIITNTQDRVIRTYELEDLLHQRG QMVEA KYKVLDMVNKA AWKNVCTDSDGLYVCGASTKAHSLYIWESNTGSLIKILHGK GEALLDV QWHPTRPIILSIAQGTVSMWTQAHVENWSAFAPEFQELEENEKYIEKENEFD MEDADE DMTSKSQDAEDDIIDVVNVRPEDYLASSDEEDCNIMKPARNLESGPLWYIPVPL DIENAD NDQRLPDDIRNLDNLLNPAWVASAKGYELSPKSK
Source	E.coli
Target Names	rbbp-5
Protein Names	Recommended name: Uncharacterized WD repeat-containing protein F21H12.1
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.