



Recombinant Putative cytosolic 5'-nucleotidase 3 (F25B5.3)

Product Code	CSB-EP600990CXY-B
Abbreviation	F25B5.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09315
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSNKVARRLG KCLFVSGRRF ESRQSILQLR TETLTDTPLS ATLDQSQFSM FKAAEIVNAA AACAEAECIE QLKKTDDVPL LMNYLLGEEQ ILVADPTAVA AKLRKMVVGG AGKTVVISDF DYTLRFANE QGERLSTTHG VFDDNVMLRK PELGQKFVDL KNKYYPPIEFS PNLTMEEKIP HMEKWWGTSH SLIVNEKFSK NTIEDFVRQS RIVFKDGAED FIALDAHNI PLVIFSAGIG NIIIEYFLQQK LGAIPRNTHF ISNMILFDED DNACAFSEPL IHTFCKNSSV IQKETSFFHD IAGRVNVILL GDSMGDIHMD VGVERDGPTL KVGYYNGSLD DTAALQHYEE VYDIVLIHDP TLNVAQKIVD IINSSH
Source	E.coli
Target Names	F25B5.3
Protein Names	Recommended name: Putative cytosolic 5'-nucleotidase 3 EC= 3.1.3.5 Alternative name(s): Putative pyrimidine 5'-nucleotidase
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.