



Recombinant Uncharacterized homeobox protein R06F6.6 (R06F6.6)

Product Code	CSB-MP601015CXY
Abbreviation	R06F6.6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09602
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MHRKGGQTLR LYTHFLLHLS PHHHSFLDSI WQRKMVHPSS AFHPYTRPIP TTSGHDLPVD DSLRLLLSAE ALIALAQLQD ASKILPSYEP VKDFSPQLLP AMTPTPIAT PSIQEQPQL QSPSAPNEKS RKRRTTFSPE QATRLEAEYI GDSYMAREKR HLLAQSLKLS ENQVKTWFQN RRAKDKRDRK SENASNHTSN SRRSSPSRKS SSDSTPTPTQ ATQFDMPTQI QTASPTTAD SAIFPPTSPE SIIQKIEQFP SNQILPNFDI LQTYLQSLSS SQIPLQFVPS TPPLFDPNML QLYHVIFEFS QSETCSVSFQ RNSASNMRSL LIFVFFATIF AVDAQYDSAD SSQSEDLPDI APESVSSEEI IEVDVNIIEGS GSGDGPIIDQ VHMFGLPGP GDIRKKRGIL EILGRK
Source	Mammalian cell
Target Names	ceh-62
Protein Names	Recommended name: Uncharacterized homeobox protein R06F6.6
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.