



# Recombinant Uncharacterized homeobox protein R06F6.6 (R06F6.6)

<b>Product Code</b>	CSB-YP601015CXY
<b>Abbreviation</b>	R06F6.6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q09602
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MHRKGGQTLR LYTHFLLHLS PHHHSFLDSI WQRKMVHPSS AFHPYTRPIP TTSGHDL PVD DSLRLLLSAE ALIALAQLQD ASKILPSYEP VKDFSPQLLP AMTPTPIAT PSIQEQPQL QSPSAPNEKS RKRRTTFSPE QATRLEAEYI GDSYMAREKR HLLAQSLKLS ENQVKTWFQN RRAKDKRDRK SENASNHTSN SRRSSPSRKS SSDSTPTPTQ ATQFDMPTQI QTASPPTTAD SAIFPPTSPE SIIQKIEQFP SNQILPNFDI LQTYLQSLSS SQIPLQFVPS TPPLFDPNML QLYHVIFEFS QSETCSVSFQ RNSASNMRSL LIFVFFATIF AVDAQYDSAD SSQSEDLPDI APESVSSEEI IEVDVNIIEGS GSGDGPIIDQ VHMFGLPGP GDIRKKRGIL EILGRK
<b>Source</b>	Yeast
<b>Target Names</b>	ceh-62
<b>Protein Names</b>	Recommended name: Uncharacterized homeobox protein R06F6.6
<b>Expression Region</b>	1-416
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.