



Recombinant Uncharacterized protein R07B1.8 (R07B1.8)

Product Code	CSB-MP601016CXY
Abbreviation	R07B1.8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09608
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTHLLLLLFF SCVQLTTGQF YCADLIQQCA KSAAEYTEQI LQYKESAFKS CVRRPACHTE RKIFDECFDA SVSATHISPP QESTSNDGET RKVTQPPMVK YNSLLKFFTE KSMNFRSALD QCFVRSFPVP KRNFFGPSIL DEDAAAYARAI YQFDLADRLW GLPELSVTRP SLDTLGVCR T QNTALRVFGS GISRIARASD PKKNNLTSSC MLDEDEITCY RQALDLNSEY IQLIYNRDYA LRACIQNLRQ QSVCRMNDKS RLRSCLCGVR EQYDNDVQAG ILQCVKSKSP HIPAMVIEMS SSLDEEPTSA KIQEQVTPAQ LTFDTPGAIV NGQCMCACEH SAKNEATRLF ANKKNNNDVE KSTQIEKKPE KQGPEIQEEV VEMETVKDEQ PPKTSAVRFK ENSPRLMQPS EAAGR VFWMN NITIT
Source	Mammalian cell
Target Names	R07B1.8
Protein Names	Recommended name: Uncharacterized protein R07B1.8
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.