



Recombinant *Schizosaccharomyces pombe* Putative glutamine amidotransferase-like protein C13C5.04 (SPAC13C5.04)

Product Code	CSB-YP601022SXV
Abbreviation	SPAC13C5.04
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09686
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDYHAILNADTPMVEITSAYGDYADMVKSLQSGSEYYSTEWNLVTKTYEVY KNPNDYP QKEDFPNINAIITGSKASATSDAPWIKKLISFVKDVLFKYPHIKIVGLCFGHQIVA KAA GVPIIQNPKGWEVSSTVVQLTENGEKFFGRKVININQMHQDMAVDVPEGFELL GSTEDCE FQIFYKPRQALTFQGHPEFSTEVVNTMVKVLRGTEVFTEQQKEEALKRSENPA DNDFLAV SIVKFLLE
Source	Yeast
Target Names	SPAC13C5.04
Protein Names	Recommended name: Putative glutamine amidotransferase-like protein C13C5.04
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.