



Recombinant Schizosaccharomyces pombe Uncharacterized protein C13G6.09 (SPAC13G6.09)

Product Code	CSB-EP601034SXV
Abbreviation	SPAC13G6.09
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09787
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MEDVDLGFLS QEKLDEKDYL DIECSRVGGA PLFIRKNDAA FLNESLENSF EFLMQLYAPK NSEISYHRIL YIFINRDGDS QTAGWTRGVK VFREQARETD EIDFQLITSK VTKLIDSSEI EVDAEPKEEN KKEIPEKLN VKVDTENPSA EPYTKAKGDV SFLKFQKRLS RAPDQIMRY Y HATSNEFPGL WCNNECIPSS IPNCACGAKR QLEFQILPTL ISSMNIDHSA KNALDWGILS IYVCSASCDL ANCGLAEEELC WLQTFSEQGI IAPE
Source	E.coli
Target Names	SPAC13G6.09
Protein Names	Recommended name: Uncharacterized protein C13G6.09
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.