



# Recombinant *Schizosaccharomyces pombe* Probable translation initiation factor eIF-2B subunit delta (tif224)

<b>Product Code</b>	CSB-MP601048SXV
<b>Abbreviation</b>	tif224
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q09924
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGFSAEQAKK DGKDQSPVSE SSSVGGTSPA TASSVSPNE PKLSGKEAKA LKKARKQASR RAKAEAAAAN NPPGVSEEKK VAIPNKNSNQ QKKASKQNPQ NSPETDANLQ EKKIFEEKQV SIFSHLDWRR RRTTENIPKD IHPAVIRLGL KLANYKIFGS NQRCIDLLKT FKIVIQDYQT PYGTTLRHL TTHINSQIAY LVSTRPLSIS MGNAIRFLKL EISVLDIDLT DDEGKELLE KIDSYIRDRI IIAGQVIVQA ATEKIQDGDV ILTYLHSSTV NDVLIHAKNV GKKFRVVVD SRPEFEGRVC LKLLTEHGIE CTYVMISALS YIMQEVTKIF LGGHAMLNSG ALYSRAGTSL ISLLGHESNV PVIACCESYK FTERIQLDSL VYNELAPGDQ LVNMGVDDFE EKPGVLNWK SVKNLKLKLSL KYDVTPPRLI TVCVCEMGLL PSTSVPAIIN EFKQVYA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	tif224
<b>Protein Names</b>	Recommended name: Probable translation initiation factor eIF-2B subunit delta Alternative name(s): eIF-2B GDP-GTP exchange factor subunit delta
<b>Expression Region</b>	1-467
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.