



Recombinant Uncharacterized protein C53C9.2 (C53C9.2)

Product Code	CSB-EP601049CXY
Abbreviation	C53C9.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09936
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGEEDWSDTQ QQPRKRWTL EQLKGGNTFLS QQAGTNKFET QKGMTAVGMP RWNITKDKKQ GYIAPDQRSE NVLRVQCGTN QYASQKGETP IGASRFQVPK VTYKKEWETI LDKEGEKIIP KQAGDYGLAS QAGEVSMGGH RNQVALIRGR LPHDRRTHGV LCFQNGTNLF ASQTGMSAPP GLGAVRQATQ KIEGLELGED ILRRGTEFTP WYSGQNKFAT QAGSGGFLKV RDVLPHTVGG KDIEEELKQK SEGIVPLQSG TNKLASQRGM TGFGTPRNTQ LLAGWKKEWI EDYEAALKEW EETKPPGSAS SVDPPFGHYKK KFEERESSRQ SEIDSQSVKA SEPVEPEPEE EEEEEEEEEKI EEPAAKEEEEE EEEEEEEEEE EELEEEEEEE EEEEEDE
Source	E.coli
Target Names	C53C9.2
Protein Names	Recommended name: Uncharacterized protein C53C9.2
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.