



# Recombinant Uncharacterized protein F12A10.8 (F12A10.8)

<b>Product Code</b>	CSB-EP601050CXY
<b>Abbreviation</b>	F12A10.8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q09949
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MHPQREVIKF RKGRFYKNPT ISSFKRIVPS NRPVAVDVNG FPQRHRRHYV QAGMFANAKS SLPSSANVSF QNSDDNLSTS RGRSASPTPI RKFSNFPRSN ISKTWDFSDI GVKKDQQTRN QLPPMKRLNS EEEEEQQGKT KTTKETIGA STVGRSVKHV LPWLKKWGSQ DAESNMKPMR IERLKVSP TL DQSINRFSSR LRRIDSIRPI PLQNHTGTTI GTVEKALVAV ALAKTKSEY QKNPVFDYLI ERSDWLLYRE EFNDKIDYLP QDPFYSHDIV SSDFEEYEQE IEEARDPYYL EIGSEEKQAP IKASASKPRF DEANFHNE DP LNLDSQEDFF ETSYLPKNAD ESKKHQKSML PSFGDYTTKS DDERQKNRLV SFDTDSSVPR TSRKRNHQEA LSPSKSNPDY SPFNYDLFSP PSVPPKPS SS SSNYSIW
<b>Source</b>	E.coli
<b>Target Names</b>	F12A10.8
<b>Protein Names</b>	Recommended name: Uncharacterized protein F12A10.8
<b>Expression Region</b>	1-437
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.