



Recombinant *Aspergillus terreus* Flap endonuclease 1 (fen1)

Product Code	CSB-YP601421APT
Abbreviation	fen1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0CBS0
Product Type	Recombinant Protein
Immunogen Species	<i>Aspergillus terreus</i> (strain NIH 2624 / FGSC A1156)
Purity	>85% (SDS-PAGE)
Sequence	MGIKQLYQVI SENAPDAIKA GDIKNHFGRK VAIDASMSIY SFLIAVRSEG QQLMSDTGET TSHLMGMFYR TLRMVDNGIK PLYVFDGAPP KLKSGELAKR TARKTEATEA HEEAKETGTA EDVEKFSRRT VRVTRHNAE CKKLLKLMGI PYIDAPTEAE AQCAVLARAG KVYAAASEDM DTLCFEAPIL LRHLTFSEQR KEPIQEIHLS RALEGLDMDR AKFIDLCILL GCDYLEPIPK VGPNTALKLI RDHGSLEKVV EYIQNDPKKK YVIPEDWPYQ DARELFLHPD VRDANDPECD FKWEAPDIEG LVEFLVKDKG FNEDRVNGA ARLQKNLKA QQSRLGFFK PVARSDDEKA TLKRKHDEKL QEQQKRKKEE AKAKKEAKAR PRGAG
Source	Yeast
Target Names	fen1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.