



Recombinant *Oryza sativa* subsp. *japonica* Putative protein phosphatase 2C 63 (Os07g0516100, LOC_Os07g33230)

Product Code	CSB-MP603538OFG
Abbreviation	Os07g0516100, LOC_Os07g33230
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0D629
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASSQQAVRE TGRGRASSSS AGGRKVTFGY HIVEGKTPHG MEDLHVAEFR RLDDGNEDGG DGGSTAVTAI LINGETLAVA NVGDSRAVAF DVRAGRAQQL SVDHEPLRER DAIEHCGGFV TEIHGDVPRV DAQLATSRAF GDRQIKEHIS SDPNVTIEDV GGRRRRWWHG ARRPRRRGV EECFLQRLAE YAICIHGASM EHRNFHVHSS HVHYFEIPPF KEV
Source	Mammalian cell
Target Names	Os07g0516100
Protein Names	Recommended name: Putative protein phosphatase 2C 63 Short name= OsPP2C63 EC= 3.1.3.16
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.