



Recombinant *Oryza sativa* subsp. *japonica* Catalase isozyme B (CATB)

| | |
|--------------------------|---|
| Product Code | CSB-YP603539OFG |
| Abbreviation | CATB |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q0D9C4 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MDPYKHRPSS GSNSTFWTTN SGAPVWNNNS ALTVGERGPI LLEDYHLIEK LAQFDRERIP ERVVHARGAS AKGFFEVTHTD ISHLTCADFL RAPGVQTPVI VRFSTVVHER GSPETLRDPR GFAVKFYTRE GNFDLVGNM PVFFIRDGMK FPDMVHAFKP SPKTNMQENW RIVDFFSHHP ESLHMFSFLF DDVGIPLNYR HMEGFGVNTY TLINKDGKPH LVKFHWKPTC GVKCLLDDEA VTGGTCHSH ATKDLTDSIA AGNYPEWKLY IQTIDPDHED RFDLFDPLDVT KTWPEDIPL QPVGRMVLNK NIDNFFAENE QLAFCPAIV PGIHYSDDKL LQTRIFSYAD TQRHRLGPNY LMLPVNAPKC AYHNNHHDGS MNFMHRDEEV NYFPSRFDAA RHAEKVPIPP RVLTGCREKC VIDKENNFQQ AGERYRSFDP ARQDRFLQRW VDALS DPRIT HELRGIWISY WSQC DASLGQ KLASRLNLKP NM |
| Source | Yeast |
| Target Names | CATB |
| Protein Names | Recommended name: Catalase isozyme B Short name= CAT-B EC= 1.11.1.6 |
| Expression Region | 1-492 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.