



Recombinant Mouse Thymosin beta-4 (Tmsb4x)

| | |
|--------------------------|---|
| Product Code | CSB-MP023934MO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P20065 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MLLPATMSDK PDMAEIEKFD KSKLKKKTETQ EKNPLPSKET IEQEKQAGES |
| Source | Mammalian cell |
| Target Names | Tmsb4x |
| Protein Names | Recommended name: Thymosin beta-4 Short name= T beta 4 Cleaved into the following chain: 1. Hematopoietic system regulatory peptide Alternative name(s): Seraspenide |
| Expression Region | 1-50 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Cytoplasmic domain |
| Target Details | This gene encodes an actin sequestering protein which plays a role in regulation of actin polymerization. The protein is also involved in cell proliferation, migration, and differentiation. This gene escapes X inactivation and has a homolog on chromosome Y. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |