



# Recombinant Bovine Synaptosomal-associated protein 29 (SNAP29)

<b>Product Code</b>	CSB-YP603793BO
<b>Abbreviation</b>	SNAP29
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0II86
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSAYPRSYNPFDEDAEDEDARPAPWSDSRDLADGPGAPADRQQALRQEVLR RAEATAAST GRSLSLMYESERIGVVSAEELVRQRGALERTEKMVDKMEQDLKTSQKHINSIK SVFGGLV NYFRSKPAETPSAQNGTLTPQPSGRLKDAINSSKEQEAQYQASHPNLRKLQD SDSIPGGA GSAVSSEAYPRNPHLRACHQRIDSNLDELSVGLGRLKDIALGIQTEIDEQDDILD RLTSK VDKLDVSITSTERKVRQL
<b>Source</b>	Yeast
<b>Target Names</b>	SNAP29
<b>Protein Names</b>	Recommended name: Synaptosomal-associated protein 29 Short name= SNAP-29 Alternative name(s): Soluble 29 kDa NSF attachment protein Vesicle-membrane fusion protein SNAP-29
<b>Expression Region</b>	1-258
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene, a member of the SNAP25 gene family, encodes a protein involved in multiple membrane trafficking steps. Two other members of this gene family, SNAP23 and SNAP25, encode proteins that bind a syntaxin protein and mediate synaptic vesicle membrane docking and fusion to the plasma membrane. This protein binds tightly to multiple syntaxins and is localized to intracellular membrane structures rather than to the plasma membrane. While the protein is mostly membrane-bound, a significant fraction of it is found free in the cytoplasm. Use of multiple polyadenylation sites has been noted for this gene.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.