



Recombinant *Oryza sativa* subsp. japonica Auxin-responsive protein IAA5 (IAA5)

Product Code	CSB-MP603820OFG
Abbreviation	IAA5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JKG7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSPPLEPHDY IGLSAAAASP TPSSSSCSSS PNPGGGARGP RLTLRLGLPG SESPEREVVA AGLTLGPLPP TTTKAASKRA FPDSSPRHGA SSGSVAAAAA CQDKAAPAAA PPAAKAQVVG WPPVRNYRKN TLAASASKGK GEDKGTAEGG PLYVKVSMDG APYLRKVDLK MYSSYEDLSM ALEKMFSCFI TGQSGLRKSS NRDRLTNGSK ADALQDQEYV LTYEDKDADW MLVGDLPWDL FTTICRKLKI MRGSDAAGIA PRSIEQSGQS R
Source	Mammalian cell
Target Names	IAA5
Protein Names	Recommended name: Auxin-responsive protein IAA5 Alternative name(s): Indoleacetic acid-induced protein 5
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.