



Recombinant Pan troglodytes NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 8 (NDUFA8)

Product Code	CSB-YP603879EQV
Abbreviation	NDUFA8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0MQB1
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	>85% (SDS-PAGE)
Sequence	PGIVELPTL EELKVDEVKI SSAVLKAAAH HYGAQCDKPN KEFMLCRWEE KDP RRCL E EG KLVNKCALDF FRQIKRHCAE PFTEYWT CID YTGQQLFRHC RKQQA K FDEC VLDKLGWVRP DLGELSKVTK VKTDRPLPEN PYHSRPRPDP SPEIEGDLQP ATHGSRFYFW TK
Source	Yeast
Target Names	NDUFA8
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 8 Alternative name(s): Complex I-19kD Short name= CI-19kD NADH-ubiquinone oxidoreductase 19 kDa subunit
Expression Region	2-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.