



Recombinant Human Uncharacterized protein C2orf80 (C2orf80)

Product Code	CSB-EP603916HU
Abbreviation	C2orf80
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0P641
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MERRLIKKEM KKLLGDYIGI RLRENEFDPK GRRQLTFLDD MAHYDLAISV ALQWLDPSED LTWLEWEELK IPLHGRPIYP NRREREMIL SSYAGILMNS IPIEEVFKIY GADSSADSGT IKVPRVSSLC LSLHPFAMLT APKAAAYARK QSVKSRKVTT NKNATSISAK EANATEWKSS QRFSDTQPKH KVT
Source	E.coli
Target Names	C2orf80
Protein Names	Recommended name: Uncharacterized protein C2orf80
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.