



Recombinant Frankia alni Lipoyl synthase (lipA)

Product Code	CSB-BP603946FAAC
Abbreviation	lipA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q0RFE8
Product Type	Recombinant Protein
Immunogen Species	Frankia alni (strain ACN14a)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MSAVAPEGRP LLRLERNAQ TPIERKPPWI RTRMRTGPEY SDVKGLVRAA GLHTVCEEAG CPNIYECWED REATFLIGGD VCTRRCDFCQ IDSGRPTPLD RDEPRRVAES VRTMGLRYAT VTGVARDDLA DGGSWLYGET VRQIHALSAG TGVEVLIPDF GGRADQLDEV FGAAPEVLAH NLETVPRIFK RIRPAFRYER SLDVLRQARA AGLVTKSNLI LGLGETAEEI HAAMRDLHAA GCELLTVTQY LRPTPRHHPV ERWVRPEEFL DWERVGAE LG FSGVMSGPLV RSSYRAGR LY QQAITARGE G HTAASDLPKS VLETSHTQ</pre>
Source	Baculovirus
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>