



Recombinant Escherichia coli O6:K15:H31 Bifunctional purine biosynthesis protein PurH (purH)

Product Code	CSB-EP604166EGY-B
Abbreviation	purH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TA59
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	<p> MQQRRPVRRA LLSVSDKAGI VEFAQALSAR GVELLSTGGT ARLLAEKGLP VTEVSDYTGF PEMMDGRVKT LHPKVHGGIL GRRGQDDAIM EEHQIQPIDM VVVNLYPFAQ TVAREGCSLE DAVENIDIGG PTMVRSAKN HRDVAIVVKS SDYDAIKEM DANEGSLTLA TRFNLAIKAF EHTAAYDSMI ANYFGSMVPA YHGESKEAAG RFPRTLNLNF IKKQDMRYGE NSHQQAIFY EENVKEASVA TATQVQGKAL SYNNIADTDA ALECVKEFAE PACVIVKHAN PCGVAVSNSI LDAYDRAYKT DPTSAFGGII AFNRELDAET AQAIISRQFV EVIIAPSASE EALKITAAKQ NVRVLCGQW GERAPGLDFK RVNGGLLVQD RDLGMVGAE LRVVTKRQPT EQELRDALFC WKVAKFVKSNI AIVYAKNNMT IGIGAGQMSR VYSAKIAGIK AADEGLEVKG SSMASDAFFP FRDGIDAAAA AGVTCVIQPG GSIRDDEVIA AADEHGIAML FTDMRHFRH </p>
Source	E.coli
Target Names	purH
Protein Names	Recommended name: Bifunctional purine biosynthesis protein PurH Including the following 2 domains: Phosphoribosylaminoimidazolecarboxamide formyltransferase EC= 2.1.2.3 Alternative name(s): AICAR transformylase IMP cyclohydrolase
Expression Region	1-529
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.