



Recombinant Escherichia coli O6:K15:H31 50S ribosomal protein L1 (rpIA)

Product Code	CSB-BP604168EGY
Abbreviation	rpIA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TA81
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MAKLTKRMRV IREKVDATKQ YDINEAIALL KELATAKFVE SVDVAVNLGI DARKSDQNV R GATVLPHTG RSVRVAVFTQ GANAEAAKAA GAELVGMEDL ADQIKKGEMN FDVVIASPDA MRVVGQLGQV LGPRGLMPNP KVGTVTPNVA EAVKNAKAGQ VRYRNDKNGI IHTTIGKVDF DADKLKENLE ALLVALKKAK PTQAKGVYIK KVSISTTMGA GVAVDQAGLS ASVN
Source	Baculovirus
Target Names	rpIA
Protein Names	Recommended name: 50S ribosomal protein L1
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.