



Recombinant Escherichia coli O6:K15:H31 Bifunctional protein GImU (glmU)

Product Code	CSB-MP604179EGY
Abbreviation	glmU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TAX9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLNNAMSVVI LAAGKGTRMY SDLPKVLHTL AGKAMVQHVI DAANELGAAH VHLVYGHGGD LLKQALKDDN LNWVLQAEQL GTGHAMQQAA PFFADDEDIL MLYGDVPLIS VETLQHLRDA KPQGGIGLLT VKLDDPTGYG RITRENGKVT GIVEHKDATD EQRRIQEINT GILIANGADM KRWLAKLTNN NAQGEYYITD IALAYQEGR EIVAVHPQRL SEVEGVNNRL QLSRLERVYQ SEQAEKLLLA GVMLRDPARF DLRGTLTHGR DVEIDTNVII EGNVTLGHRV KIGTGCVIKN SVIGDDCEIS PYTVVEDANL AAACITGPFA RLRPGAELLE GAHVGNFVEM KKARLGKGSK AGHLTYLGDA EIGDNVNIGA GTITCNYDGA NKFKTIIGDD VFGSDTQLV APVTVGKGAT IAAGTTVTRN VGENALAIR VPQTQKEGWR RPVKKK
Source	Mammalian cell
Target Names	glmU
Protein Names	Recommended name: Bifunctional protein GImU Including the following 2 domains: UDP-N-acetylglucosamine pyrophosphorylase EC= 2.7.7.23 Alternative name(s): N-acetylglucosamine-1-phosphate uridyltransferase Glucosamine-1-phosphate
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.