



Recombinant Escherichia coli O6:K15:H31 5'/3'-nucleotidase SurE (surE)

Product Code	CSB-BP604205EGY
Abbreviation	surE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TEB4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRILLSNDDG VHAPGIQTLA KALREFADVQ VVAPDRNRSG ASNSLTLESS LRTFTFENGDI AVQMGTPTD CVYLGVNALM RRPDIVVSG INAGPNLGDD VIYSGTVAAA MEGRHLGFPA LAVSLDGHKH YDTAAAVTCS ILRALCKEPL RTGRILNINV PDLPLDQIKG IRVTRCGTRH PADQVIPQQD PRGNTLYWIG PPGGKCDAGP GTDFAAVDEG YVSITPLHVD LTAHSAQDVV SDWLNSVGVG TQW
Source	Baculovirus
Target Names	surE
Protein Names	Recommended name: 5'/3'-nucleotidase SurE EC= 3.1.3.5 EC= 3.1.3.6 Alternative name(s): Exopolyphosphatase EC= 3.6.1.11 Nucleoside monophosphate phosphohydrolase
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.