



Recombinant Escherichia coli O6:K15:H31 Anaerobic nitric oxide reductase flavorubredoxin (norV)

Product Code	CSB-MP604206EGY
Abbreviation	norV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TEH0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MSIVVKNNIH WVGQRDWEVR DFHGTEYKTL RGSSYNSYLI REEKNVLIDT VDHKFSREFV QNLRNEIDLA DIDYIVINHA EEDHAGALTE LMAQIPDTPI YCTANAIDSI NGHHHHPewn FNVVKTGDTL DIGNGKQLIF VETPMLHWPD SMMTYLTGDA VLFSNDAFGQ HYCDEHLFND EVDQTELFEQ CQRYANILT PFSRLVTPKI TEILGFNLPV DMIATSHGVV WRDNPTQIVE LYLKWAADYQ EDRITIVYDT MSNNTRMMAD AIAQGIAETD PRVAVKIFNV ARSDKNEILT NVFRSKGVLV GTSTMNNVMM PKIAGLVEEM TGLRFRNKRA SAFGSHGWSG GAVDRLSTRL QDAGFEMSLS LKAKWRPDQD ALELCREHGR EIARQWALAP LPQSTVNTVV EEETSAATTA DLGPRMQCSV CQWIYDPAKG EPMQDVAPGT PWSEVPDNFL CPECSLGKDV FDELASEAK
Source	Mammalian cell
Target Names	norV
Protein Names	Recommended name: Anaerobic nitric oxide reductase flavorubredoxin Short name= FIRd Short name= FlavoRb
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.