



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At2g30100, chloroplastic (At2g30100)

Product Code	CSB-BP604494DOA
Abbreviation	At2g30100
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0WNN7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	EMGLSRSVEL DQFITSEEEE GEAEIEIGEGF FEAIEELERM TREPSDILEE MNHRLSSREL QLMLVYFAQE GRDSWCTLEV FEWLKKNRV DEEIMELMVS IMCGWVKKLI EDECNAHQVF DLLIEMDCVG LKPGFSMMDK VIALYCEMGK KESAVLFVKE VLRRRDGFY SVVGGGGSEG RKGGPVGYLA WKFMVDGDYR KAVDMVMELR LSGLKPEAYS YLIAMTAIVK ELNSLGKTLR ELKRFARAGF VAEIDDHDRV LIEKYQSETL SRGLQLATWA VEEGQENDSI IGVVHERLLA MYICAGRGPE AEKQLWKMKL AGREPEADLH DIVMAICASQ KEVNAVSRLL TRVEFMGSQR KKKTLWLLR GYVKGGHFEE AAETLVSMID SGLHPEYIDR VAVMQGMTRK IQRPRDVEAY MSLCKRLFDA GLVGPCLVYM YIDKYKLWIV KMM
Source	Baculovirus
Target Names	At2g30100
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At2g30100, chloroplastic
Expression Region	51-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.