



Recombinant Arabidopsis thaliana AT-rich interactive domain-containing protein 5 (ARID5)

Product Code	CSB-MP604496DOA
Abbreviation	ARID5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0WNR6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MMADTEMQEQ DVPSGTKGVV EEPSELEKDL NSIERPKSPV PEDTTHTLDS DVHLSDAPIA NQIEANEEVG GQNSVDGRNG DVDQSEKKIT SDGGQEETTL GESNPLKGDV SSPHVPEESV KKWKTWLLSD AEAREVDEAG APQDQEAFIK EVEAFNKENF LEFKAPKFYG QPLNCLKLWR AVIKLGGYDV VTTSKLWRQV GESFHPPKTC TTVSWTFRIF YEKALLEYEK HLRQNGELNL PGSASLPSSG IEKEASSHQA SGSGRTRRDA AARAMQGWHS QRLLGSGEVT EPIVKEKGLN STPKQKNLKN IGVQKQKTTT GMDLVFSHES EKQSTAEVID VGPPADWVKI NVRETKDCFE IFALVPGLLR EEVRVQSDPA GRLVIAGQPE QLDNPWGITP FKKVVNFPAR IDPLHTSAVV SLHGRLFVRV PFEQ
Source	Mammalian cell
Target Names	ARID5
Protein Names	Recommended name: AT-rich interactive domain-containing protein 5 Short name= ARID domain-containing protein 5
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.