



Recombinant Schizosaccharomyces pombe UPF0012 hydrolase C26A3.11 (SPAC26A3.11)

Product Code	CSB-MP604564SXV
Abbreviation	SPAC26A3.11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10166
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNSKFFGLVQ KGTRSFPSL NFCYTRNIMS VSASSLVPKD FRAFRIGLVQ LANTKDKSEN LQLARLKVLE AAKNGSNVIV LPEIFNSPYG TGYFNQYAEP IEESSPSYQA LSSMAKDTKT YLFGGSIPER KDGKLYNTAM VFDPSGKLIA VHRKIHLFDI DIPGGVSFRE SDSLSPGDAM TMVDTEYGKF GLGICYDIRF PELAMIAARN GCSVMIYPGA FNLSTGPLHW ELLARARAVD NEMFVACCAP ARDMNADYHS WGHSTVDPF GKVIATTEK PSIVYADIDP SVMSTARNSV PIYTQRRFDV YSEVLPALKK EE
Source	Mammalian cell
Target Names	SPAC26A3.11
Protein Names	Recommended name: UPF0012 hydrolase C26A3.11 EC= 3.5.-.-
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.