



Recombinant *Schizosaccharomyces pombe* Meiotically up-regulated gene 152 protein (mug152)

Product Code	CSB-EP604581SXV-B
Abbreviation	mug152
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10274
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MRSMKPPGFS SDLMDEHSVD LLNGSILAAE NPSKREVAQD VPGFERKPTK VRKPRVKWTE KETNDLLRGC QIHGVGNWKK ILLDERFHFT NRSPNDLKDR FRTILPEDYK KFYPNAKTHM GRPQKIPHTV GLSKSTRKER KQFTPEEDER LLEGFFLHGP CWTRISKDAN LGLQNRSTD LRDRFRNAFP ERYAAAGFKL KNNPGNRSKY YQNNMVNDAT TPNDSSSTEA AAAVAVAVA VAASNPNASP QQTTEQPASD ELLDWP HHNL PSQFFTSQRN PNYSTDSFLL GQSLSDPFNH TLQSFHPYES LFSAGQPPSL PISPSTSQNS VQPFPSIQQ PPLHLEPPLS SNTLNSSTLP QPNSTDFNTF PPLPSTPRIS SEDIPWDNRG
Source	E.coli
Target Names	teb1
Protein Names	Recommended name: Meiotically up-regulated gene 152 protein
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.