



Recombinant Dual specificity mitogen-activated protein kinase kinase mek-2 (mek-2)

Product Code	CSB-EP604623CXY-B
Abbreviation	mek-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10664
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSSGKRRNPL GLSLPPTVNE QSESGEATAE EATATVPLEE QLKKLGLTEP QTQRLSEFLQ VKEGIKELSE DMLQTEGELG HGNGGVVVKC VHRKTGVIMA RKL VHLEIKP SVRQQIVKEL AVLHKCNSPF IVGFYGA FVD NNDISICMEY MDGLSLDIVL KKVGR LPEKF VGRISVAVVR GLTYLKDEIK ILHRDVKPSN MLVNSNGEIK LCDFGVSGML IDSMANSFVG TRSYMAPERL TGSHYTISSD IWSFGLSLVE LLIGRYPVPA PSQAEYATMF NVAENEIELA DSLEEPNYHP PSNPASMAIF EMLDYIVNGP PPTLPKRFFT DEVIGFVSKC LRKLP SERAT LKSLTADVFF TQYADHDDQG EFAV FVKGTI NLPKLNP
Source	E.coli
Target Names	mek-2
Protein Names	Recommended name: Dual specificity mitogen-activated protein kinase kinase mek-2 Short name= MAP kinase kinase mek-2 EC= 2.7.12.2
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.